



DCAF8 rabbit pAb

Catalog No	YP-Ab-09593
Isotype	IgG
Reactivity	Human;Mouse;Rat
Applications	WB
Gene Name	DCAF8 H326 WDR42A
Immunogen	Synthesized peptide derived from human DCAF8 AA range: 65-115
Specificity	This antibody detects endogenous levels of human DCAF8
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Polyclonal, Rabbit,IgG
Purification	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Dilution	WB 1:1000-2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	DCAF8;DDB1- and CUL4-associated factor 8;H326;WD repeat-containing protein 42A;DCAF8 H326 WDR42A;DDB1- and CUL4-associated factor 8 (WD repeat-containing protein 42A)
Observed Band	67-73 kDa
Calculated Molecular Weight	597 aa, 67 kDa
Cell Pathway	Nucleus . Cytoplasm . It shuttles between the nucleus and the cytoplasm. Nuclear import is mediated by KPNA1 and KPNB1 under the regulation of nuclear GTPase RAN. Nuclear export to the cytoplasm is XPO1 dependent. .
Tissue Specificity	
Function	
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images