



SETD4 rabbit pAb

Catalog No	YP-Ab-09594
Isotype	IgG
Reactivity	Human;Mouse
Applications	WB
Gene Name	SETD4 C21orf18 C21orf27
Immunogen	Synthesized peptide derived from human SETD4 AA range: 193-243
Specificity	This antibody detects endogenous levels of human SETD4
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Polyclonal, Rabbit,IgG
Purification	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Dilution	WB 1:1000-2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	C21orf18;C21orf27;SET domain containing 4;SETD 4;SETD4;SETD4 C21orf18 C21orf27;SET domain-containing protein 4
Observed Band	40 kDa
Calculated Molecular Weight	440 aa, 50 kDa
Cell Pathway	Cytoplasm, cytosol . Nucleus .
Tissue Specificity	
Function	
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images