



ACSF4 rabbit pAb

Catalog No	YP-Ab-09595
Isotype	IgG
Reactivity	Human;Mouse
Applications	WB
Gene Name	AASDH ACSF4 U26 HSPC318
Immunogen	Synthesized peptide derived from human ACSF4 AA range: 673-723
Specificity	This antibody detects endogenous levels of human ACSF4
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Polyclonal, Rabbit,IgG
Purification	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Dilution	WB 1:1000-2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	ACSF4;AASDH ACSF4 U26 HSPC318;Acyl-CoA synthetase family member 4 (EC 6.2.1.-) (Amino adipate-semialdehyde dehydrogenase) (Protein NRPS998)
Observed Band	
Calculated Molecular Weight	
Cell Pathway	
Tissue Specificity	Ubiquitously expressed in adult tissues.
Function	
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images

