



MRGBP rabbit pAb

Catalog No	YP-Ab-09598
Isotype	IgG
Reactivity	Human;Mouse
Applications	WB
Gene Name	MRGBP C20orf20
Immunogen	Synthesized peptide derived from human MRGBP AA range: 99-149
Specificity	This antibody detects endogenous levels of human MRGBP
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Polyclonal, Rabbit,IgG
Purification	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Dilution	WB 1:1000-2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	Eaf7;MRG/MORF4L-binding protein;MRG15BP;MRG-binding protein;MRGBP;MRGBP C20orf20
Observed Band	30 kDa
Calculated Molecular Weight	204 aa, 22 kDa
Cell Pathway	Nucleus.
Tissue Specificity	
Function	
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images
