



ARL14 rabbit pAb

Catalog No	YP-Ab-09599
Isotype	IgG
Reactivity	Human;Mouse
Applications	WB
Gene Name	ARL14 ARF7
Immunogen	Synthesized peptide derived from human ARL14 AA range: 67-117
Specificity	This antibody detects endogenous levels of human ARL14
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Polyclonal, Rabbit,IgG
Purification	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Dilution	WB 1:1000-2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	ADP ribosylation factor like 14;ADP-ribosylation factor 7;ADP-ribosylation factor-like protein 14;ARF7;ARL14;ARL14_HUMAN;FLJ22595;OTTHUMP213879;ARL14 ARF7;ADP-ribosylation factor-like protein 14 (ADP-ribosylation factor 7)
Observed Band	
Calculated Molecular Weight	22kD
Cell Pathway	Cytoplasmic vesicle . Colocalizes with MHC II-containing cytoplasmic vesicles.
Tissue Specificity	Expressed in immature dendritic cells.
Function	
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images