



ZN644 rabbit pAb

Catalog No	YP-Ab-09600
Isotype	IgG
Reactivity	Human;Rat;Mouse;
Applications	WB
Gene Name	ZNF644 KIAA1221 ZEP2
Immunogen	Synthesized peptide derived from human ZN644 AA range: 674-724
Specificity	This antibody detects endogenous levels of human ZN644
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Polyclonal, Rabbit,IgG
Purification	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Dilution	WB 1:1000-2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	ZN644;ZNF644 KIAA1221 ZEP2;Zinc finger protein 644 (Zinc finger motif enhancer-binding protein 2) (Zep-2)
Observed Band	
Calculated Molecular Weight	
Cell Pathway	Nucleus.
Tissue Specificity	Expressed in liver, placenta, retina and retinal pigment epithelium.
Function	
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images
