



FXL16 rabbit pAb

Catalog No	YP-Ab-09602
Isotype	IgG
Reactivity	Human;Mouse;Rat
Applications	WB
Gene Name	FBXL16 C16orf22 FBL16
Immunogen	Synthesized peptide derived from human FXL16 AA range: 395-445
Specificity	This antibody detects endogenous levels of human FXL16
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Polyclonal, Rabbit,IgG
Purification	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Dilution	WB 1:1000-2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	FXL16;FBXL16 C16orf22 FBL16;F-box/LRR-repeat protein 16 (F-box and leucine-rich repeat protein 16)
Observed Band	
Calculated Molecular Weight	
Cell Pathway	
Tissue Specificity	
Function	
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images