



TMM95 rabbit pAb

Catalog No	YP-Ab-09606
Isotype	IgG
Reactivity	Human;Mouse
Applications	WB
Gene Name	TMEM95 UNQ9390/PRO34281
Immunogen	Synthesized peptide derived from human TMM95 AA range: 103-153
Specificity	This antibody detects endogenous levels of human TMM95
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Polyclonal, Rabbit,IgG
Purification	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Dilution	WB 1:1000-2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	TMM95;TMEM95 UNQ9390/PRO34281;Transmembrane protein 95
Observed Band	
Calculated Molecular Weight	
Cell Pathway	Cytoplasmic vesicle, secretory vesicle, acrosome membrane ; Single-pass type I membrane protein . Following the acrosome reaction, relocalizes to the equatorial segment. .
Tissue Specificity	
Function	
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images