



TM249 rabbit pAb

Catalog No	YP-Ab-09607
Isotype	IgG
Reactivity	Human;Rat;Mouse;
Applications	WB
Gene Name	TMEM249 C8orfK29
Immunogen	Synthesized peptide derived from human TM249 AA range: 133-183
Specificity	This antibody detects endogenous levels of human TM249
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Polyclonal, Rabbit,IgG
Purification	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Dilution	WB 1:1000-2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	TM249;TMEM249 C8orfK29;Transmembrane protein 249
Observed Band	
Calculated Molecular Weight	
Cell Pathway	Cell projection, cilium, flagellum membrane ; Multi-pass membrane protein . Predominantly located in the principal piece of the sperm tail. .
Tissue Specificity	
Function	
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images
