



ZFAN5 rabbit pAb

Catalog No	YP-Ab-09608
Isotype	IgG
Reactivity	Human;Mouse;Rat
Applications	WB
Gene Name	ZFAND5 ZA20D2 ZNF216
Immunogen	Synthesized peptide derived from human ZFAN5 AA range: 84-134
Specificity	This antibody detects endogenous levels of human ZFAN5
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Polyclonal, Rabbit,IgG
Purification	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Dilution	WB 1:1000-2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	AN1-type zinc finger protein 5;ZA2D2;ZFAN5_HUMAN;Zfand5;ZFAND5A;zinc finger A2 domain containing 2;Zinc finger A2 domain containing protein 2;Zinc finger A2 domain-containing protein 2;zinc finger AN1 type domain 5;Zinc finger protein 216;ZNF216;ZFAN5;ZFAND5 ZA20D2 ZNF216;AN1-type zinc finger protein 5 (Zinc finger A20 domain-containing protein 2) (Zinc finger protein 216)
Observed Band	
Calculated Molecular Weight	23kD
Cell Pathway	Cytoplasm .
Tissue Specificity	Highly expressed in skeletal muscle. Expressed in fetal cochlea. Also expressed in infant brain, fetal heart, pancreatic islet, melanocyte, pineal gland, placenta, corneal stroma, and parathyroid tumor. Weakly expressed or undetectable in adult brain, heart, colon, thymus, spleen, kidney, liver, small intestine, placenta, lung and peripheral blood leukocytes. Expressed in rhabdomyosarcoma RD cells (at protein level).
Function	
Background	

**matters needing attention**

Avoid repeated freezing and thawing!

Usage suggestions

This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.

Products Images