



SCC4 rabbit pAb

Catalog No	YP-Ab-09609
Isotype	IgG
Reactivity	Human;Mouse
Applications	WB
Gene Name	MAU2 KIAA0892 SCC4
Immunogen	Synthesized peptide derived from human SCC4 AA range: 558-608
Specificity	This antibody detects endogenous levels of human SCC4
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Polyclonal, Rabbit,IgG
Purification	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Dilution	WB 1:1000-2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	KIAA0892;mau 2;MAU2L;Protein MAU 2;SCC4;MAU2 KIAA0892 SCC4;MAU2 chromatid cohesion factor homolog (MAU-2) (Cohesin loading complex subunit SCC4 homolog)
Observed Band	66-69 kDa
Calculated Molecular Weight	613 aa, 69 kDa
Cell Pathway	Nucleus, nucleoplasm . Nucleus . Chromosome . Binds to chromatin from the end of mitosis until prophase.
Tissue Specificity	
Function	
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images
