



SE6L2 rabbit pAb

Catalog No	YP-Ab-09610
Isotype	IgG
Reactivity	Human;Mouse
Applications	WB
Gene Name	SEZ6L2 PSK UNQ1903/PRO4349
Immunogen	Synthesized peptide derived from human SE6L2 AA range: 106-156
Specificity	This antibody detects endogenous levels of human SE6L2
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Polyclonal, Rabbit,IgG
Purification	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Dilution	WB 1:1000-2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	SE6L2;SEZ6L2 PSK UNQ1903/PRO4349;Seizure 6-like protein 2
Observed Band	
Calculated Molecular Weight	
Cell Pathway	Cell membrane ; Single-pass type I membrane protein . Endoplasmic reticulum membrane ; Single-pass type I membrane protein . Detected on cell surface of lung-cancer.
Tissue Specificity	
Function	
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images
