



INT12 rabbit pAb

Catalog No	YP-Ab-09612
Isotype	IgG
Reactivity	Human;Mouse;Rat
Applications	WB
Gene Name	INTS12 PHF22 SBBI22
Immunogen	Synthesized peptide derived from human INT12 AA range: 398-448
Specificity	This antibody detects endogenous levels of human INT12
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Polyclonal, Rabbit,IgG
Purification	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Dilution	WB 1:1000-2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	Int12;PHF22;SBBI22;INTS12 PHF22 SBBI22;Integrator complex subunit 12 (Int12) (PHD finger protein 22)
Observed Band	49 kDa
Calculated Molecular Weight	462 aa, 49 kDa
Cell Pathway	Nucleus .
Tissue Specificity	
Function	
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images
