



FA92B rabbit pAb

Catalog No	YP-Ab-09614
Isotype	IgG
Reactivity	Human;Mouse
Applications	WB
Gene Name	FAM92B
Immunogen	Synthesized peptide derived from human FA92B AA range: 187-237
Specificity	This antibody detects endogenous levels of human FA92B
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Polyclonal, Rabbit,IgG
Purification	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Dilution	WB 1:1000-2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	FA92B;FAM92B;Protein FAM92B
Observed Band	
Calculated Molecular Weight	
Cell Pathway	Cytoplasm, cytoskeleton, microtubule organizing center, centrosome, centriole . Cytoplasm, cytoskeleton, cilium basal body . Extensive colocalization with CBY1 at mother centrioles. .
Tissue Specificity	Restricted to certain tissues, most prominently expressed in multiciliated tissues.
Function	
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images
