



RN216 rabbit pAb

Catalog No	YP-Ab-09616
Isotype	IgG
Reactivity	Human;Mouse
Applications	WB
Gene Name	RNF216 TRIAD3 UBCE7IP1 ZIN
Immunogen	Synthesized peptide derived from human RN216 AA range: 196-246
Specificity	This antibody detects endogenous levels of human RN216
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Polyclonal, Rabbit,IgG
Purification	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Dilution	WB 1:1000-2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	RN216;RNF216 TRIAD3 UBCE7IP1 ZIN;E3 ubiquitin-protein ligase RNF216 (EC 6.3.2.-) (RING finger protein 216) (Triad domain-containing protein 3) (Ubiquitin-conjugating enzyme 7-interacting protein 1) (Zinc finger protein inhibiting NF-
Observed Band	
Calculated Molecular Weight	
Cell Pathway	Cytoplasm.
Tissue Specificity	Ubiquitous, with the highest levels of expression in testis and peripheral blood leukocytes.
Function	
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images
