



PKHH1 rabbit pAb

Catalog No	YP-Ab-09617
Isotype	IgG
Reactivity	Human;Mouse
Applications	WB
Gene Name	PLEKHH1 KIAA1200
Immunogen	Synthesized peptide derived from human PKHH1 AA range: 335-385
Specificity	This antibody detects endogenous levels of human PKHH1
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Polyclonal, Rabbit,IgG
Purification	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Dilution	WB 1:1000-2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	PKHH1;PLEKHH1 KIAA1200;Pleckstrin homology domain-containing family H member 1 (PH domain-containing family H member 1)
Observed Band	
Calculated Molecular Weight	
Cell Pathway	
Tissue Specificity	
Function	
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images
