



S35G5 rabbit pAb

Catalog No	YP-Ab-09619
Isotype	IgG
Reactivity	Human;Rat;Mouse;
Applications	WB
Gene Name	SLC35G5 AMAC AMAC1L2
Immunogen	Synthesized peptide derived from human S35G5 AA range: 47-97
Specificity	This antibody detects endogenous levels of human S35G5
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Polyclonal, Rabbit,IgG
Purification	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Dilution	WB 1:1000-2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	S35G5;SLC35G5 AMAC AMAC1L2;Solute carrier family 35 member G5 (Acyl-malonyl-condensing enzyme 1-like protein 2)
Observed Band	
Calculated Molecular Weight	
Cell Pathway	Membrane ; Multi-pass membrane protein .
Tissue Specificity	Expressed in placenta and testis.
Function	
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images
