



RPB1B rabbit pAb

Catalog No	YP-Ab-09620
Isotype	IgG
Reactivity	Human;Rat;Mouse;
Applications	WB
Gene Name	POLR2J2
Immunogen	Synthesized peptide derived from human RPB1B AA range: 18-68
Specificity	This antibody detects endogenous levels of human RPB1B
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Polyclonal, Rabbit,IgG
Purification	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Dilution	WB 1:1000-2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	RPB1B;POLR2J2;DNA-directed RNA polymerase II subunit RPB11-b1 (RNA polymerase II subunit B11-b1) (RPB11b1) (DNA-directed RNA polymerase II subunit J2)
Observed Band	
Calculated Molecular Weight	
Cell Pathway	Nucleus .
Tissue Specificity	Ubiquitously expressed.
Function	
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images
