



UBP2 rabbit pAb

Catalog No	YP-Ab-09622
Isotype	IgG
Reactivity	Human;Mouse;Rat
Applications	WB
Gene Name	USP2 UBP41
Immunogen	Synthesized peptide derived from human UBP2 AA range: 452-502
Specificity	This antibody detects endogenous levels of human UBP2
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Polyclonal, Rabbit,IgG
Purification	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Dilution	WB 1:1000-2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	Deubiquitinating enzyme 15;Deubiquitinating enzyme;EC 3 1 2 15;KIAA529; MGC131982;MGC149838;MGC74854;Ubiquitin Carboxy terminal Hydrolase 15; Ubiquitin carboxyl-terminal hydrolase 15;Ubiquitin specific peptidase 15; Ubiquitin specific processing protease 15;Ubiquitin Specific Protease 15; Ubiquitin thioesterase 15;Ubiquitin thiolesterase 15;Ubiquitin-specific-processing protease 15;UBP 15;UBP15;UBP15_HUMAN;Unph 2;UNPH 4;Unph-2;Unph2; Unph4;USP 15;Usp15;UBP2;USP2 UBP41;Ubiquitin carboxyl-terminal hydrolase 2 (EC 3.4.19.12) (41 kDa ubiquitin-specific protease) (Deubiquitinating enzyme 2) (Ubiquitin thioesterase 2) (Ubiquitin-specific-processing protease 2)
Observed Band	
Calculated Molecular Weight	67kD
Cell Pathway	Cytoplasm . Cytoplasm, perinuclear region . Localizes in the spermatid head in late-elongating spermatids in the thin area between the outer acrosomal membrane and the plasma membrane. . ; [Isoform 4]: Nucleus . Membrane ; Peripheral membrane protein . Cytoplasm . Predominantly expressed at membranes. .
Tissue Specificity	Expressed in mesangial cells of the kidney and in different types of glomerulonephritides (at protein level).

Function**Background****matters needing attention**

Avoid repeated freezing and thawing!

Usage suggestions

This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.

Products Images