



ZWINT rabbit pAb

Catalog No	YP-Ab-09623
Isotype	IgG
Reactivity	Human;Mouse;Rat
Applications	WB
Gene Name	ZWINT
Immunogen	Synthesized peptide derived from human ZWINT AA range: 165-215
Specificity	This antibody detects endogenous levels of human ZWINT
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Polyclonal, Rabbit,IgG
Purification	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Dilution	WB 1:1000-2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	HZwint 1;KNTC2AP;ZW10-interacting protein 1;Zwint 1;ZWINT1;ZWINT;ZW10 interactor (ZW10-interacting protein 1) (Zwint-1)
Observed Band	35 kDa
Calculated Molecular Weight	277 aa, 31 kDa
Cell Pathway	Nucleus. Chromosome, centromere, kinetochore. Localizes to kinetochores from late prophase to anaphase.
Tissue Specificity	
Function	
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images
