



FNDC4 rabbit pAb

Catalog No	YP-Ab-09625
Isotype	IgG
Reactivity	Human;Mouse
Applications	WB
Gene Name	FNDC4 FRCP1 UNQ6389/PRO21134
Immunogen	Synthesized peptide derived from human FNDC4 AA range: 33-83
Specificity	This antibody detects endogenous levels of human FNDC4
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Polyclonal, Rabbit,IgG
Purification	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Dilution	WB 1:1000-2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	FNDC4;FNDC4 FRCP1 UNQ6389/PRO21134;Fibronectin type III domain-containing protein 4 (Fibronectin type III repeat-containing protein 1)
Observed Band	
Calculated Molecular Weight	
Cell Pathway	Membrane ; Single-pass type I membrane protein . Secreted . The N-terminus is probably cleaved to release a secreted extracellular portion of the protein. .
Tissue Specificity	Up-regulated in colon cells of patients with inflammatory bowel disease (IBD).
Function	
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images
