



C1QT7 rabbit pAb

Catalog No	YP-Ab-09626
Isotype	IgG
Reactivity	Human;Mouse
Applications	WB
Gene Name	C1QTNF7 CTRP7
Immunogen	Synthesized peptide derived from human C1QT7 AA range: 53-103
Specificity	This antibody detects endogenous levels of human C1QT7
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Polyclonal, Rabbit,IgG
Purification	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Dilution	WB 1:1000-2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	C1q and tumor necrosis factor related protein 7;C1QT7 HUMAN;C1QTNF7; Complement C1q tumor necrosis factor related protein 7;Complement C1q tumor necrosis factor-related protein 7;CTRP 7;ZACRP7;C1QT7;C1QTNF7 CTRP7
Observed Band	
Calculated Molecular Weight	31kD
Cell Pathway	Secreted .
Tissue Specificity	
Function	
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images
