



GBP6 rabbit pAb

Catalog No	YP-Ab-09627
Isotype	IgG
Reactivity	Human;Rat;Mouse;
Applications	WB
Gene Name	GBP6
Immunogen	Synthesized peptide derived from human GBP6 AA range: 243-293
Specificity	This antibody detects endogenous levels of human GBP6
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Polyclonal, Rabbit,IgG
Purification	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Dilution	WB 1:1000-2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	Guanylate-binding protein 6;Guanine nucleotide-binding protein 6;GTP-binding protein 6;GBP-6;EC:3.6.5.-;GBP6;Guanylate-binding protein 6 (GTP-binding protein 6) (GBP-6) (Guanine nucleotide-binding protein 6)
Observed Band	70 kDa
Calculated Molecular Weight	633 aa, 72 kDa
Cell Pathway	Cytoplasmic vesicle .
Tissue Specificity	
Function	
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images
