



CAN14 rabbit pAb

Catalog No	YP-Ab-09630
Isotype	IgG
Reactivity	Human;Rat;Mouse;
Applications	WB
Gene Name	CAPN14
Immunogen	Synthesized peptide derived from human CAN14 AA range: 136-186
Specificity	This antibody detects endogenous levels of human CAN14
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Polyclonal, Rabbit,IgG
Purification	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Dilution	WB 1:1000-2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	CAN14;CAPN14;Calpain-14 (EC 3.4.22.-) (Calcium-activated neutral proteinase 14) (CANP 14)
Observed Band	
Calculated Molecular Weight	
Cell Pathway	
Tissue Specificity	Not expressed in tissues tested.
Function	
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images
