



ECM2 rabbit pAb

Catalog No	YP-Ab-09632
Isotype	IgG
Reactivity	Human;Mouse
Applications	WB
Gene Name	ECM2
Immunogen	Synthesized peptide derived from human ECM2 AA range: 101-151
Specificity	This antibody detects endogenous levels of human ECM2
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Polyclonal, Rabbit,IgG
Purification	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Dilution	WB 1:1000-2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	Matrix glycoprotein SC1/ECM2;ECM2;Extracellular matrix protein 2 (Matrix glycoprotein SC1/ECM2)
Observed Band	
Calculated Molecular Weight	699 aa, 8 kDa
Cell Pathway	Secreted, extracellular space, extracellular matrix .
Tissue Specificity	Expressed predominantly in adipose tissue as well as female-specific organs such as mammary gland, ovary, and uterus.
Function	
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images
