



GNPTG rabbit pAb

Catalog No	YP-Ab-09635
Isotype	IgG
Reactivity	Human;Mouse
Applications	WB
Gene Name	GNPTG C16orf27 GNPTAG CAB56184 LP2537
Immunogen	Synthesized peptide derived from human GNPTG AA range: 216-266
Specificity	This antibody detects endogenous levels of human GNPTG
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Polyclonal, Rabbit,IgG
Purification	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Dilution	WB 1:1000-2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	C16orf27;CAB56184;GlcNAc-1-phosphotransferase subunit gamma;GNPTAG; LP2537;GNPTG;GNPTG C16orf27 GNPTAG CAB56184 LP2537;N-acetylglucosamine-1-phosphotransferase subunit gamma (GlcNAc-1-phosphotransferase subunit gamma) (UDP-N-acetylglucosamine-1-phosphotransferase subunit gamma)
Observed Band	34 kDa, 68~70 kDa
Calculated Molecular Weight	34 kDa
Cell Pathway	Secreted . Golgi apparatus .
Tissue Specificity	Widely expressed.
Function	
Background	
matters needing attention	Avoid repeated freezing and thawing!

**Usage suggestions**

This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images
