



REC16 rabbit pAb

Catalog No	YP-Ab-09638
Isotype	IgG
Reactivity	Human;Rat;Mouse;
Applications	WB
Gene Name	
Immunogen	Synthesized peptide derived from human REC16 AA range: 8-58
Specificity	This antibody detects endogenous levels of human REC16
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Polyclonal, Rabbit,IgG
Purification	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Dilution	WB 1:1000-2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	REC16;HERV-K_10p14 provirus Rec protein
Observed Band	
Calculated Molecular Weight	
Cell Pathway	Cytoplasm . Nucleus, nucleolus . Shuttles between the nucleus and the cytoplasm. When in the nucleus, resides in the nucleolus (By similarity). .
Tissue Specificity	
Function	
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images
