



ACADS rabbit pAb

Catalog No	YP-Ab-09640
Isotype	IgG
Reactivity	Human;Mouse;Rat
Applications	WB
Gene Name	ACADS
Immunogen	Synthesized peptide derived from human ACADS AA range: 347-397
Specificity	This antibody detects endogenous levels of human ACADS
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Polyclonal, Rabbit,IgG
Purification	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Dilution	WB 1:1000-2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	ACAD3;ACADS;Butyryl CoA dehydrogenase;SCAD;Short-chain specific acyl-CoA dehydrogenase;mitochondrial (SCAD) (EC 1.3.8.1) (Butyryl-CoA dehydrogenase)
Observed Band	42 kDa
Calculated Molecular Weight	412 aa, 44 kDa
Cell Pathway	Mitochondrion matrix .
Tissue Specificity	
Function	
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images

