



ROGDI rabbit pAb

Catalog No	YP-Ab-09643
Isotype	IgG
Reactivity	Human;Mouse;Rat
Applications	WB
Gene Name	ROGDI
Immunogen	Synthesized peptide derived from human ROGDI AA range: 148-198
Specificity	This antibody detects endogenous levels of human ROGDI
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Polyclonal, Rabbit,IgG
Purification	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Dilution	WB 1:1000-2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	Protein rogdi homolog;ROGDI;rogdi homolog (Drosophila)
Observed Band	32 kDa
Calculated Molecular Weight	287 aa, 32 kDa
Cell Pathway	Nucleus envelope . Cell junction, synapse, presynapse . Cell projection, axon . Perikaryon . Cell projection, dendrite . Cytoplasmic vesicle, secretory vesicle, synaptic vesicle . Detected primarily at presynaptic sites on axons, and to a lesser degree in soma and dendrites. Not detected at post-synaptic sites. .
Tissue Specificity	Widely expressed with highest levels in spinal cord, brain, heart and bone marrow. Also expressed in fetal brain and liver.
Function	
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images