



FA59B rabbit pAb

Catalog No	YP-Ab-09645
Isotype	IgG
Reactivity	Human;Mouse
Applications	WB
Gene Name	FAM59B KIAA2038 HRIHFB2063
Immunogen	Synthesized peptide derived from human FA59B AA range: 385-435
Specificity	This antibody detects endogenous levels of human FA59B
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Polyclonal, Rabbit,IgG
Purification	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Dilution	WB 1:1000-2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	FA59B_HUMAN;FAM 59B;FAM59B;Family with sequence similarity 59 member B;FLJ375;HRIHFB263;KIAA238;Protein FAM59B;FA59B;FAM59B KIAA2038 HRIHFB2063
Observed Band	
Calculated Molecular Weight	93kD
Cell Pathway	
Tissue Specificity	
Function	
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images

