



CNPY4 rabbit pAb

Catalog No	YP-Ab-09649
Isotype	IgG
Reactivity	Human;Mouse
Applications	WB
Gene Name	CNPY4 PSEC0237 UNQ1909/PRO4354
Immunogen	Synthesized peptide derived from human CNPY4 AA range: 167-217
Specificity	This antibody detects endogenous levels of human CNPY4
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Polyclonal, Rabbit,IgG
Purification	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Dilution	WB 1:1000-2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	CNPY4;CNPY4 PSEC0237 UNQ1909/PRO4354;Protein canopy homolog 4
Observed Band	
Calculated Molecular Weight	
Cell Pathway	Secreted .
Tissue Specificity	
Function	
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images
