



TDRD1 rabbit pAb

Catalog No	YP-Ab-09652
Isotype	IgG
Reactivity	Human;Mouse
Applications	WB
Gene Name	TDRD1
Immunogen	Synthesized peptide derived from human TDRD1 AA range: 803-853
Specificity	This antibody detects endogenous levels of human TDRD1
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Polyclonal, Rabbit,IgG
Purification	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Dilution	WB 1:1000-2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	hPCL3;MGC131698;MGC149712;MGC149713;MGC23929;MTF2L1;OTTHUMP21996;OTTHUMP21998;PCL3;PHD finger protein 19;Phf19;PHF19 protein;PHF19_HUMAN;Polycomb like 3;Polycomb like protein 3;Polycomb-like protein 3;TDRD19B;Tudor domain containing 19B;TDRD1;Tudor domain-containing protein 1 (Cancer/testis antigen 41.1) (CT41.1)
Observed Band	130kD
Calculated Molecular Weight	66 kDa; 66kD
Cell Pathway	Cytoplasm . Component of the meiotic nuage, also named P granule, a germ-cell-specific organelle required to repress transposon activity during meiosis. Also present in chromatoid body (By similarity). .
Tissue Specificity	Testis and ovary specific. Also expressed in several cancers.
Function	
Background	
matters needing attention	Avoid repeated freezing and thawing!

**Usage suggestions**

This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images
