



SRXN1 rabbit pAb

Catalog No	YP-Ab-09654
Isotype	IgG
Reactivity	Human;Mouse
Applications	WB
Gene Name	SRXN1 C20orf139 SRX
Immunogen	Synthesized peptide derived from human SRXN1 AA range: 40-90
Specificity	This antibody detects endogenous levels of human SRXN1
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Polyclonal, Rabbit,IgG
Purification	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Dilution	WB 1:1000-2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	C20orf139;EC:1.8.98.2;SRX;SRX1;Sulfiredoxin 1;SRXN1;SRXN1 C20orf139 SRX;Sulfiredoxin-1 (EC 1.8.98.2)
Observed Band	14 kDa
Calculated Molecular Weight	14 kDa
Cell Pathway	Cytoplasm .
Tissue Specificity	Widely expressed with highest levels in kidney, lung, spleen and thymus.
Function	
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images
