



FRMD3 rabbit pAb

Catalog No	YP-Ab-09657
Isotype	IgG
Reactivity	Human;Mouse
Applications	WB
Gene Name	FRMD3 EPB41L4O
Immunogen	Synthesized peptide derived from human FRMD3 AA range: 238-288
Specificity	This antibody detects endogenous levels of human FRMD3
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Polyclonal, Rabbit,IgG
Purification	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Dilution	WB 1:1000-2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	FRMD3;FRMD3 EPB41L4O;FERM domain-containing protein 3 (Band 4.1-like protein 4O) (Ovary type protein 4.1) (4.1O)
Observed Band	
Calculated Molecular Weight	
Cell Pathway	Membrane ; Single-pass membrane protein .
Tissue Specificity	Ovary-specific.
Function	
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images
