



PBX1 rabbit pAb

Catalog No	YP-Ab-09660
Isotype	IgG
Reactivity	Human;Mouse
Applications	WB
Gene Name	PBX1 PRL
Immunogen	Synthesized peptide derived from human PBX1 AA range: 162-212
Specificity	This antibody detects endogenous levels of human PBX1
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Polyclonal, Rabbit,IgG
Purification	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Dilution	WB 1:1000-2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	Homeobox protein PBX1;Homeobox protein PRL;pre B cell leukemia homeobox 1;Pre-B-cell leukemia transcription factor 1;PRL;PBX1;PBX1 PRL;Pre-B-cell leukemia transcription factor 1 (Homeobox protein PBX1) (Homeobox protein PRL)
Calculated Molecular Weight	430 aa, 47 kDa
Cell Pathway	Nucleus .
Tissue Specificity	Expressed in the kidney. Expressed in the endothelial cells of the glomeruli and interstitium (at protein level) (PubMed:28270404). Expressed in all tissues except in cells of the B and T lineage. Expressed strongly in kidney and brain (PubMed:28270404).
Function	
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images