



MGN2 rabbit pAb

Catalog No	YP-Ab-09661
Isotype	IgG
Reactivity	Human;Mouse
Applications	WB
Gene Name	MAGOHB MAGOH2
Immunogen	Synthesized peptide derived from human MGN2 AA range: 55-105
Specificity	This antibody detects endogenous levels of human MGN2
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Polyclonal, Rabbit,IgG
Purification	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Dilution	WB 1:1000-2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	High mobility group (nonhistone chromosomal) protein 17;high mobility group nucleosomal binding domain 2;High mobility group nucleosome-binding domain-containing protein 2;High mobility group protein N2;HMG17;HMGN2; HMGN2_HUMAN;MGC5629;Non histone chromosomal protein HMG 17;Non-histone chromosomal protein HMG-17;Nonhistone chromosomal protein HMG 17;nonhistone chromosomal protein hmg-17;MGN2;MAGOHB MAGOH2;Protein mago nashi homolog 2
Observed Band	
Calculated Molecular Weight	9 KD; 9kD
Cell Pathway	Nucleus .
Tissue Specificity	Ubiquitous.
Function	
Background	
matters needing attention	Avoid repeated freezing and thawing!

**Usage suggestions**

This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images
