



SRTD4 rabbit pAb

Catalog No	YP-Ab-09662
Isotype	IgG
Reactivity	Human;Mouse
Applications	WB
Gene Name	SERTAD4
Immunogen	Synthesized peptide derived from human SRTD4 AA range: 203-253
Specificity	This antibody detects endogenous levels of human SRTD4
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Polyclonal, Rabbit,IgG
Purification	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Dilution	WB 1:1000-2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	ATD4;IFT139;JBTS11;Nbla1696;NPHP12;Putative protein product of Nbla1696; SRTD4;Tetratricopeptide repeat protein 21B;THM1;TPR repeat protein 21B; TT21B_HUMAN;Ttc21b;SERTAD4;SERTA domain-containing protein 4
Observed Band	
Calculated Molecular Weight	151kD
Cell Pathway	
Tissue Specificity	
Function	
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images
