



F125A rabbit pAb

Catalog No	YP-Ab-09663
Isotype	IgG
Reactivity	Human;Mouse;Rat
Applications	WB
Gene Name	FAM125A CFBP
Immunogen	Synthesized peptide derived from human F125A AA range: 11-61
Specificity	This antibody detects endogenous levels of human F125A
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Polyclonal, Rabbit,IgG
Purification	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Dilution	WB 1:1000-2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	F125A;FAM125A CFBP;Multivesicular body subunit 12A (CIN85/CD2AP family-binding protein) (ESCRT-I complex subunit MVB12A) (Protein FAM125A)
Observed Band	
Calculated Molecular Weight	
Cell Pathway	Cytoplasm. Nucleus. Endosome. Cytoplasm, cytoskeleton, microtubule organizing center, centrosome. Late endosome membrane ; Peripheral membrane protein . Colocalizes with F-actin. Some fraction may be nuclear.
Tissue Specificity	Ubiquitously expressed except in skeletal muscle.
Function	
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images
