



GTPB6 rabbit pAb

Catalog No	YP-Ab-09666
Isotype	IgG
Reactivity	Human;Mouse
Applications	WB
Gene Name	GTPBP6 PGPL
Immunogen	Synthesized peptide derived from human GTPB6 AA range: 410-460
Specificity	This antibody detects endogenous levels of human GTPB6
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Polyclonal, Rabbit,IgG
Purification	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Dilution	WB 1:1000-2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	GTPB6;GTPBP6 PGPL;Putative GTP-binding protein 6 (Pseudoautosomal GTP-binding protein-like)
Observed Band	
Calculated Molecular Weight	
Cell Pathway	
Tissue Specificity	Ubiquitously expressed.
Function	
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images
