



EMC4 rabbit pAb

Catalog No	YP-Ab-09668
Isotype	IgG
Reactivity	Human;Mouse
Applications	WB
Gene Name	EMC4 TMEM85 HSPC184 PIG17
Immunogen	Synthesized peptide derived from human EMC4 AA range: 91-141
Specificity	This antibody detects endogenous levels of human EMC4
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Polyclonal, Rabbit,IgG
Purification	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Dilution	WB 1:1000-2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	EMC4;ER membrane protein complex subunit 4;PIG17;TMEM85; transmembrane protein 85;EMC4 TMEM85 HSPC184 PIG17;ER membrane protein complex subunit 4 (Cell proliferation-inducing gene 17 protein) (Transmembrane protein 85)
Observed Band	20 kDa
Calculated Molecular Weight	
Cell Pathway	Endoplasmic reticulum membrane ; Multi-pass membrane protein . Could also be a single-pass transmembrane protein with cytosolic N-terminus and luminal C-terminus. .
Tissue Specificity	Isoform 1 is expressed in brain and heart. Isoform 2 is expressed in heart.
Function	
Background	
matters needing attention	Avoid repeated freezing and thawing!

**Usage suggestions**

This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images
