



PK1L2 rabbit pAb

Catalog No	YP-Ab-09669
Isotype	IgG
Reactivity	Human;Mouse
Applications	WB
Gene Name	PKD1L2 KIAA1879 PC1L2
Immunogen	Synthesized peptide derived from human PK1L2 AA range: 1505-1555
Specificity	This antibody detects endogenous levels of human PK1L2
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Polyclonal, Rabbit,IgG
Purification	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Dilution	WB 1:1000-2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	PK1L2;PKD1L2 KIAA1879 PC1L2;Polycystic kidney disease protein 1-like 2 (PC1-like 2 protein) (Polycystin-1L2)
Observed Band	
Calculated Molecular Weight	
Cell Pathway	Membrane ; Multi-pass membrane protein .
Tissue Specificity	Expressed in all tissues tested including brain, placenta, mammary gland, testis, lung and liver. Highest expression in skeletal muscle. Isoform 2 is expressed in heart and kidney.
Function	
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images
