



P5I13 rabbit pAb

Catalog No	YP-Ab-09675
Isotype	IgG
Reactivity	Human;Mouse;Rat
Applications	WB
Gene Name	TP53I13 DSCP1
Immunogen	Synthesized peptide derived from human P5I13 AA range: 36-86
Specificity	This antibody detects endogenous levels of human P5I13
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Polyclonal, Rabbit,IgG
Purification	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Dilution	WB 1:1000-2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	P5I13;TP53I13 DSCP1;Tumor protein p53-inducible protein 13 (Damage-stimulated cytoplasmic protein 1)
Observed Band	
Calculated Molecular Weight	
Cell Pathway	Cell membrane ; Single-pass type I membrane protein ; Extracellular side . Cytoplasm . Associates with unknown subcellular structures in the cytoplasm.
Tissue Specificity	Expressed in heart, placenta, skeletal muscle, testis, brain and lung.
Function	
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images
