



BARX1 rabbit pAb

Catalog No	YP-Ab-09677
Isotype	IgG
Reactivity	Human;Mouse
Applications	WB
Gene Name	BARX1
Immunogen	Synthesized peptide derived from human BARX1 AA range: 159-209
Specificity	This antibody detects endogenous levels of human BARX1
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Polyclonal, Rabbit,IgG
Purification	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Dilution	WB 1:1000-2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	BarH like homeobox 1;BarX 1;BARX homeobox 1;Homeobox protein BarH like 1; BARX1;Homeobox protein BarH-like 1
Observed Band	
Calculated Molecular Weight	27kD; 27kD(
Cell Pathway	Nucleus .
Tissue Specificity	Widely expressed. Expressed at higher levels in testis and heart. Detected in craniofacial tissue and adult iris, but not in lymphocytes, fibroblasts, choroid retina, retinal pigment epithelium, kidney, or fetal liver.
Function	
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images
