



MYRIP rabbit pAb

Catalog No	YP-Ab-09680
Isotype	IgG
Reactivity	Human;Mouse;Rat
Applications	WB
Gene Name	MYRIP SLAC2C
Immunogen	Synthesized peptide derived from human MYRIP AA range: 591-641
Specificity	This antibody detects endogenous levels of human MYRIP
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Polyclonal, Rabbit,IgG
Purification	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Dilution	WB 1:1000-2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	MYRIP;MYRIP SLAC2C;Rab effector MyRIP (Exophilin-8) (Myosin-VIIa- and Rab-interacting protein) (Synaptotagmin-like protein lacking C2 domains C) (SlAC2-c) (Slp homolog lacking C2 domains c)
Observed Band	
Calculated Molecular Weight	
Cell Pathway	Cytoplasm . Cytoplasm, perinuclear region . Cytoplasmic vesicle, secretory vesicle . In presynaptic and postsynaptic areas in photoreceptor cells and in the basal microvilli of retinal pigment epithelium cells. Associated with melanosomes. Colocalizes with actin filaments. .
Tissue Specificity	Detected in brain, skin, heart, adrenal medulla, pancreas, intestine, liver, kidney, muscle and testis.
Function	
Background	
matters needing attention	Avoid repeated freezing and thawing!

**Usage suggestions**

This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.

Products Images