



NTNG2 rabbit pAb

Catalog No	YP-Ab-09682
Isotype	IgG
Reactivity	Human;Mouse
Applications	WB
Gene Name	NTNG2 KIAA1857 LMNT2 UNQ9381/PRO34206
Immunogen	Synthesized peptide derived from human NTNG2 AA range: 384-434
Specificity	This antibody detects endogenous levels of human NTNG2
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Polyclonal, Rabbit,IgG
Purification	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Dilution	WB 1:1000-2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	KIAA1857;Laminet-2;LMNT2;Netrin G2;Netrin-G2;NTNG2;NTNG2 KIAA1857 LMNT2 UNQ9381/PRO34206;Netrin-G2 (Laminet-2)
Observed Band	50 kDa
Calculated Molecular Weight	60 kDa
Cell Pathway	Cell membrane ; Lipid-anchor, GPI-anchor ; Extracellular side .
Tissue Specificity	
Function	
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images
