



TOPZ1 rabbit pAb

Catalog No	YP-Ab-09684
Isotype	IgG
Reactivity	Human;Mouse;Rat
Applications	WB
Gene Name	TOPAZ1 C3orf77
Immunogen	Synthesized peptide derived from human TOPZ1 AA range: 583-633
Specificity	This antibody detects endogenous levels of human TOPZ1
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Polyclonal, Rabbit,IgG
Purification	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Dilution	WB 1:1000-2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	TOPZ1;TOPAZ1 C3orf77;Testis- and ovary-specific PAZ domain-containing protein 1
Observed Band	
Calculated Molecular Weight	
Cell Pathway	Cytoplasm, cytosol .
Tissue Specificity	
Function	
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images
