



POK7 rabbit pAb

Catalog No	YP-Ab-09685
Isotype	IgG
Reactivity	Human;Rat;Mouse;
Applications	WB
Gene Name	
Immunogen	Synthesized peptide derived from human POK7 AA range: 450-500
Specificity	This antibody detects endogenous levels of human POK7
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Polyclonal, Rabbit,IgG
Purification	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Dilution	WB 1:1000-2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	DKFZp547O146;Factor binding IST protein 1;Factor that binds to inducer of short transcripts protein 1;FBI-1;FBI1;HIV-1 1st-binding protein 1;HIV-1 inducer of short transcripts binding protein;HIV-1 inducer of short transcripts-binding factor 1;Leukemia/lymphoma-related factor;LRF;lymphoma related factor; MGC99631;POK erythroid myeloid ontogenic factor;Pokemon;POZ and Krueppel erythroid myeloid ontogenic factor;TIP21;TTF-I-interacting peptide 21; ZBT7A_HUMAN;ZBTB7;ZBTB7A;Zinc finger and BTB domain containing 7A; Zinc finger and BTB domain-containing protein 7A;Zinc finger protein 857A;Zinc finger- and BTB domain-containing protein 7;ZNF857A;POK7;HERV-K_1q23.3 provirus ancestral Pol protein (HERV-K(C1a) Pol protein) (HERV-K110 Pol protein) (HERV-K18 Pol protein) [Includes: Reverse transcriptase (EC 2.7.7.49); Ribonuclease H (RNase H) (EC 3.1.2
Observed Band	
Calculated Molecular Weight	61kD; 61kD
Cell Pathway	
Tissue Specificity	
Function	



Background

matters needing attention

Avoid repeated freezing and thawing!

Usage suggestions

This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images
