



MFS11 rabbit pAb

Catalog No	YP-Ab-09687
Isotype	IgG
Reactivity	Human;Mouse
Applications	WB
Gene Name	MFSD11 ET
Immunogen	Synthesized peptide derived from human MFS11 AA range: 323-373
Specificity	This antibody detects endogenous levels of human MFS11
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Polyclonal, Rabbit,IgG
Purification	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Dilution	WB 1:1000-2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	MFS11;MFSD11 ET;UNC93-like protein MFSD11 (Major facilitator superfamily domain-containing protein 11) (Protein ET)
Observed Band	
Calculated Molecular Weight	
Cell Pathway	Membrane ; Multi-pass membrane protein .
Tissue Specificity	
Function	
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images
