



KTN1 rabbit pAb

Catalog No	YP-Ab-09694
Isotype	IgG
Reactivity	Human;Mouse
Applications	WB
Gene Name	KTN1 CG1 KIAA0004
Immunogen	Synthesized peptide derived from human KTN1 AA range: 485-535
Specificity	This antibody detects endogenous levels of human KTN1
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Polyclonal, Rabbit,IgG
Purification	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Dilution	WB 1:1000-2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	KTN 1;KIAA0004;CG-1 antigen;KTN1;KTN1 CG1 KIAA0004;Kinectin (CG-1 antigen) (Kinesin receptor)
Observed Band	156 kDa
Calculated Molecular Weight	1357 aa, 156 kDa
Cell Pathway	Endoplasmic reticulum membrane; Single-pass type II membrane protein. Vesicle membrane protein anchored to the endoplasmic reticulum.
Tissue Specificity	High levels in peripheral blood lymphocytes, testis and ovary, lower levels in spleen, thymus, prostate, small intestine and colon.
Function	
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images

