



OMP rabbit pAb

Catalog No	YP-Ab-09695
Isotype	IgG
Reactivity	Human;Mouse;Rat
Applications	WB
Gene Name	OMP
Immunogen	Synthesized peptide derived from human OMP AA range: 31-81
Specificity	This antibody detects endogenous levels of human OMP
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Polyclonal, Rabbit,IgG
Purification	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Dilution	WB 1:1000-2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	Olfactory marker protein;Olfactory neuronal specific protein;OMP;Olfactory marker protein (Olfactory neuronal-specific protein)
Observed Band	
Calculated Molecular Weight	18kDa; 19kD
Cell Pathway	Cytoplasm.
Tissue Specificity	Uniquely associated with mature olfactory receptor neurons.
Function	
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images
